

Ultrasonic Cleaner



DG-CL65



DG-CL100

Products Features:

Dual power switch, Digital control(1-30Min),

Stainless steel 304, Clean and safe

20 degree bevel design body touch control

Stainless steel cleaning basket and valve

Avoid oxidation of cleaning items

Ultrasonic cleaning, emulsification, mixing, extraction, dispersion, degassing, etc.

DATA Sheets

Model	DG-CL65	DG-CL100
Volume	6.5L	10L
Ultrasonic power	90W/180W Dual power switching	120W/240W Dual power switchin
Ultrasonic frequency	40KHz	
Heating power	200Wsumnred smtoihnang	200Wwumredomuohnaong
Time control	Digital control(1-30Min)	
Temperature control	Room Temperature-80°C	
Chamber volume	300×150×150mm	300×240×150mm
Chamber material	Stainless steel 304	
Outer dimension	330×180×310mm	330×270×310mm
Packing size(L×W×H)	400×255×355mm	430×345×380mm
Control panel	20 degree bevel design body touch control	
Degassing function	Avoid oxidation of cleaning items	
Water valve	Stainless steel water valve	
Cleaning basket	Stainless steel cleaning basket	
Voltage	AC100-120V,AC220V-240V	
Weight	N.W.5.5kg,G.W.6kg	N.W.7.5kg,G.W.8.3kg