



Parameter: SpO2, HR



Big LED screen with two direction display




8s auto shut down without operation.



Superior SPO2 technology offer accurate test result.



Anti-drop, waterproof and dustproof design, longer service.



Environmental silicone layer, comfortable to wear.

Parameter:

Measurement Specification

SpO2 displayed range	0%~100%
SpO2 accuracy	70%~100%±2% 0%~69%: not specified
PR displayed range	18bpm~250bpm
PR accuracy	25~250bpm±3bpm 18~24bpm, Not specified
Data update period	<30s
Data average time	8s
SpO2 sensor peak wavelength range	Red light: 660±3nm Infrared light: 940±5nm
Maximum photic output power	<30mW

