

Water Distiller

Products Features:

The filter bamboo charcoal bag

Simple operation, safe and reliable

Super large water outflow capacity

High quality stainless inner tank and water faucet.

Double temperature control protection modeContinuous
productionEasy to clean



Technical Data

Rated voltage: 220VAC 50Hz 110VAC 60Hz

Rated power: 750V Volume: 4L

Distillation volume: 1L/HNet

weight: 4.kg Gross weight: 4.8kg

Cavity size: ϕ 180x200mm

Dimensions :280x240x390mm

Packing size :290x260x460mm

NOTE: With automatic over temperature protectionfunction, When
the temperature is higher than 115degrees, the power supply will
be cut off automatically